


INSTRUCTIONS FOR USE

Product Name: DNA Selection Beads

Catalog # DSB01-D291

- Step 1.** Leave the selection beads at room temp for 30 min prior to use to equilibrate.
- Step 2.** Vortex or pipette up and down 10 times to mix the beads thoroughly.
- Step 3.** Add the first round of selection beads to the sample and mix thoroughly. Incubate at room temp for 5 min.
- Step 4.** Add the second round of selection beads to the sample and mix thoroughly. Incubate at room temp for 5 min.
- Step 5.** Spin down the tube and place the tube on a magnetic stand.
- Step 6.** When the solution is clear (5 min wait), aspirate the supernatant and discard.
- Step 7.** Add 200 μ L of 80% Ethanol while on the magnetic stand. Incubate at room temp for 30 sec. Aspirate the ethanol and discard.
- Step 8.** Repeat the previous step one more time.
- Step 9.** Remove residual ethanol with 10 μ L pipette tips and lightly air dry the selection beads (~5min).
 Over drying the beads may result in a lower recovery DNA target. Leave the tube open only part way.
- Step 10.** Remove the tube from the magnetic stand and add $\geq 20\mu$ L ddH₂O and mix thoroughly. Incubate at room temp for 5 min.



Recommended Conditions for DNA Size Selection:

Length of DNA Fragment	250-30 bp	320-420 bp	450-550 bp	550-70 bp	700-900 bp	800-1000 bp
Ratio of Beads: DNA for 1 st Round	0.80x	0.70x	0.60x	0.55x	0.50x	0.45x
Ratio of Beads: DNA for 2 nd Round	0.20x	0.20x	0.20x	0.15x	0.15x	0.15x

"x" refers to the volume of the sample DNA.

If the insert length is 250bp and the sample DNA volume is 100 μ L, the volume of beads used in the 1st round is 0.80x100 μ L=80 μ L.