

INSTRUCTIONS FOR USE

Product Name: T7 RNA Polymerase, GMP-Grade (250 U/μL)

Catalog # T702-E51B

In vitro transcription without cap-analog:

Step 1. In a centrifuge tube add the following reaction components:

2 μL (1x concentration)	10x Transcription Buffer
2 μL each (10mM each)	CTP / GTP/ ATP/ UTP (100 mM each)
5 μL	T7 RNA Polymerase (50 U/μL)
0.04 μL	Pyrophosphatase, Inorganic (1 U/μL)
0.5 μL	RNase inhibitor (40 U/μL)
Up to 18 μL	RNase free H ₂ O
2 μL (100ng-1μg)	DNA template

Note: The DNA template should be added last. The 10x Transcription Buffer contains high concentration of spermidine, which may cause precipitation of DNA templates.

- It is advisable to keep the water and buffer at room temperature before use. The reaction mix should also be prepared at room temperature because spermidine will cause precipitation of high-concentration DNA templates at low temperature.
- For transcripts less than 100 nt: the templates should be increased to 2μg.
- To ensure effective transcription of a specific region, it is recommended to cut the DNA template downstream region into blunt ends or 5' protruding ends.

Step 2. Incubate the reaction at 37°C for 2-4 h (For transcripts less than 100 nt: increase the incubating time to 4-8 h).

Step 3. To degrade the DNA template: Add 2U DNase I and incubate at 37°C for 15-30 min.

Step 4. Purification of transcripts: RNA Cleaner magnetic beads can be used to remove proteins, salt ions and other impurities. Phenol/chloroform purification method can also be used.

In vitro transcription with cap1-analog:

Step 1. In a centrifuge tube add the following reaction components:

2 μL (1x concentration)	10x Transcription Buffer
2 μL each (10mM each)	CTP / GTP/ ATP/ UTP (100 mM each)
2 μL each (10mM each)	Cap1-analog (100 mM each)
5 μL	T7 RNA Polymerase (50 U/μL)
0.04 μL	Pyrophosphatase, Inorganic (1 U/μL)

0.5 µL	RNase inhibitor (40 U/µL)
Up to 18 µL	RNase free H ₂ O
2 µL (100ng-1µg)	DNA template

Step 2. Incubate the reaction at 37°C for 2-4 h.

Step 3. To degrade the DNA template: Add 2U DNase I and incubate at 37°C for 15-30 min.

Step 4. Purification of transcripts: RNA Cleaner magnetic beads can be used to remove proteins, salt ions and other impurities. Phenol/chloroform purification method can also be used.