

INSTRUCTIONS FOR USE

Product Name: T7 High-Yield RNA Synthesis Kit

Catalog # T702-K598

DNA Template preparation: Linearized plasmids with T7 promoters or PCR products can serve as templates. Resuspend in TE buffer or RNase-free H₂O. T7 Promoter Sequence: TAATACGACTCACTATAG*GG (*G = transcription start site).

A) Plasmid Template:

Step 1. Insert the target DNA into a plasmid vector containing the T7 promoter.

Step 2. Linearize completely using restriction enzymes and purify.

- Circular plasmids yield heterogeneous RNA; linearization ensures specific-length transcripts.
- Use restriction enzymes that cut downstream of the insert (5' overhangs or blunt ends preferred).
- Purify linearized plasmids to remove enzymes/salts that may inhibit transcription.

B) PCR Product Template:

Step 1. Add the T7 promoter sequence to the 5' end of the upstream primer.

Step 2. Amplify using high-fidelity polymerase.

- Confirm PCR product specificity/concentration via electrophoresis. Use 2-5 µL per 20 µL reaction.
- Gel purification improves RNA yield and quality.

In Vitro Transcription:

Step 1. Centrifuge T7 RNA Polymerase Mix briefly; keep on ice.

Step 2. Thaw 10× Transcription Buffer and NTPs; vortex and spin down.

Step 3. Setup reaction (20 µL System):

2 µL	10× Transcription Buffer
2 µL each (10 mM each)	ATP/CTP/GTP/UTP (100 mM each)
1 µg	DNA Template
2 µL	T7 RNA Polymerase Mix
Up to 20 µL	RNase-free H ₂ O

- Assemble at room temperature (spermidine in buffer may precipitate DNA at low temps).
- For >1000nt transcripts: Prefer linearized plasmids.
- Use a PCR machine with an open lid to prevent evaporation.
- White precipitate (magnesium pyrophosphate) is harmless; add EDTA if needed.

Step 4. Mix, spin down, and incubate at 37°C for 2 hours (extend to 4-8 hours for short <100nt transcripts).

Step 5. (Optional) Add 2 µL RNase-free DNase I and incubate at 37°C for 15 min to remove DNA template.

RNA purification:

A) RNA Cleaner Magnetic Beads:

- Step 1.** Equilibrate RNA Clean Beads to room temperature (~30 min).
- Step 2.** Dilute transcription product to 50 µL with RNase-free H₂O.
- Step 3.** Add 100 µL beads (2× volume), pipette mix 6×.
- Step 4.** Incubate 5 min at RT for RNA binding.
- Step 5.** Place on magnetic stand for 5 min; discard supernatant.
- Step 6.** Wash twice with 200 µL freshly prepared 80% ethanol (30 sec per wash).
- Step 7.** Air-dry beads 5 min (avoid cracking).
- Step 8.** Resuspend in 22 µL RNase-free H₂O; incubate 5 min.
- Step 9.** Transfer 20 µL supernatant to a new tube (avoid beads).
- Step 10.** Store at -80°C if not used immediately.

B) Phenol/Chloroform Extraction:

- Step 1.** Dilute 20 µL reaction with 115 µL H₂O + 15 µL 3M NaOAc (pH 5.2).
- Step 2.** Extract with 1× phenol/chloroform (1:1) → 2× chloroform.
- Step 3.** Add 2× volume ethanol; incubate ≥30 min at -20°C.
- Step 4.** Centrifuge 15 min (max speed, 4°C); wash pellet with 500 µL 70% ethanol.
- Step 5.** Resuspend in 20 µL RNase-free H₂O; store at -80°C.

C) Lithium Chloride Precipitation (for RNA length >300 nt and ≥100 ng/µL):

- Step 1.** Add 30 µL H₂O + 30 µL 7.5M LiCl to 20 µL reaction.
- Step 2.** Incubate ≥30 min at -20°C; centrifuge as above.
- Step 3.** Wash pellet with 70% ethanol; resuspend in 20 µL H₂O.

D) Column Purification:

Dilute product to 100 µL with RNase-free H₂O; follow column protocol.

RNA Quantification:

UV Absorption: Measure A260/A280.

Dye-Based (RiboGreen): Accurate for purified/unpurified RNA (unaffected by free nucleotides).

RNA Quality Assessment:

Gel Electrophoresis: Use agarose or polyacrylamide gels to check size/integrity.

Agilent 2100 Bioanalyzer: Requires minimal RNA; high-quality samples show sharp peaks.

Notes:

- Use RNase-free tips/tubes and avoid RNase introduction during handling.