

INSTRUCTIONS FOR USE

Product Name: mRNA Vaccinia Capping Enzyme, GMP-Grade (10 U/ μ L)**Catalog # VAC01-E59B**

Cap1 capping reaction (20 μ L reaction system):

This step is suitable for capping reaction of 10 μ g RNA (≥ 100 nt) and can be amplified according to experimental needs.

Step 1. In a 1.5 ml centrifuge tube add 10 μ g of RNA and dilute to 9.5 μ L with nuclease-free water.

Step 2. Heat for 5 minutes at 65°C.

Step 3. Place the centrifuge tube on ice for 5 minutes.

Step 4. Add the following components in sequence:

9.5 μ L	Denatured RNA
2.0 μ L	10x Capping Buffer
0.5 μ L	Murine RNase inhibitor (40U U/ μ L)
1.0 μ L	GTP (10 mM)
1.0 μ L	SAM (10 mM, fresh)
5.0 μ L	Vaccinia Capping Enzyme (10 U/ μ L)
1.0 μ L	Cap 2' -O-Methyltransferase (50 U/ μ L)

Note: 10x Capping Buffer: 0.5 M Tris-HCl, 50 mM KCl, 10 mM MgCl₂, 10 mM DTT pH 8.0 at 25°C.

Step 5. Incubate at 37°C for 2 hours.

Step 6. RNA capping is completed, next experiments can be performed.