


SPECIFICATION SHEET

Product Name: Der f 26**Catalog # DF26-A35H**

Product Description: Recombinant *Dermatophagoides farinae* allergen Der f 26 was expressed in *E. Coli* and purified by affinity chromatography using a His tag.

Storage and Stability: Store at -20°C. To avoid repeated handling and multiple freeze/thaw cycles, aliquot product into smaller quantities.

Alternative Names	: Myosin light chain (MLC)
Source	: <i>Dermatophagoides farinae</i>
Accession Numbers	: GenBank: KM009996
Host	: <i>E. Coli</i>
Formulation	: Protein stored in 1x PBS
Molecular Weight	: Calculated: 20 kDa Observed: ~20 kDa
Purity	: The purity was determined to be $\geq 90\%$ by densitometry
Gel Image	



References

1. An S, et al. *Dermatophagoides farinae* allergens diversity identification by proteomics. Mol Cell Proteomics. 2013 Jul;12(7):1818-28.
2. Suojalehto H, et al. Occupational Mite Allergy and Asthma: An EAACI Task Force Report. Allergy. 2025 Sep;80(9):2484-2500.

Version #: 1.0

To place your order, please contact us by phone 1-778-326-0223 or 1-888-606-3424 (Toll free)
or by email: orders@signalchemdx.com or info@signalchemdx.com - www.signalchemdx.com