


SPECIFICATION SHEET

Product Name: Ory c 3
Catalog # OC03-A25H

Product Description: Recombinant *Oryctolagus cuniculus* allergen Ory c 3 was expressed in *P. pastoris* and purified by affinity chromatography using a His tag.

Storage and Stability: Store at -20°C. To avoid repeated handling and multiple freeze/thaw cycles aliquot product into smaller quantities.

Alternative Names : Rabbit lipophilin, rabbit secretoglobulin
Source : *Oryctolagus cuniculus* (rabbit)
Accession Numbers : Genbank: [AAG42806](#)
 UniProt: [Q9GK63](#)
Host : *P. pastoris*
Formulation : Protein stored in 20 mM Tris, 300 mM NaCl, 10% Glycerol, pH 7.5
Molecular Weight : Calculated: 18.3 kDa
 Observed: ~18.0 kDa
Purity : The purity was determined to be ≥ 90% by densitometry

Gel Image

References

1. Hilger C, et al. Identification and isolation of a Fel d 1-like molecule as a major rabbit allergen. *J Allergy Clin Immunol.* 2014 Mar;133(3):759-66.
2. Kler S, et al. Evaluation and characterization of Ory c 3, a major rabbit allergen. *Clinical & Translational Allergy (CTA)* 2018, 8(Suppl 1):P30

Version: 1.0